

Liquid Reagents – ready to use

# LDH-P (Lactate Dehydrogenase - P) IFCC

2 Reagents

**Diagnostic reagent for quantitative in vitro determination of lactate dehydrogenase (LDH) in human serum or plasma on photometric systems.**

Ref.No.	Kit Size	Content
DIA010341	62.5 ml	2x25 ml R1 + 1x12.5 ml R2
DIA010342	5 x 20 ml	4x20 ml R1 + 1x20 ml R2
DIA010343	5 x 25 ml	4x25 ml R1 + 1x25 ml R2

Additionally offered:

DIA040012	1 x 3 mL	Diacal Auto (Calibrator)
DIA030012	1 x 5 mL	Diacon N (Control Normal)
DIA030022	1 x 5 mL	Diacon P (Control Abnormal)

### TEST PARAMETERS

<b>Method:</b>	UV, Kinetic, Decreasing Reaction Optimized DGKC
<b>Wavelength:</b>	Hg 334 nm, Hg 365 nm, 340 nm
<b>Temperature:</b>	25°C, 30°C, 37°C
<b>Sample:</b>	Serum, heparin- or EDTA-plasma,
<b>Linearity:</b>	up to 3059 U/L
<b>Sensitivity:</b>	The lower limit of detection is 5 U/L

### REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
<b>Reagent 1:</b>	
Pyruvate	0.60 mmol/L
Phosphate	50 mmol/L
<b>Reagent 2:</b>	
NADH	0.18 mmol/L
Good's buffer, pH 9.6	

### REAGENT PREPARATION

#### Substrate Start:

Reagents are ready for use.

#### Sample Start:

Mix 4 parts of Reagent 1 with 1 part of Reagent 2 (= Working Reagent).

### REAGENT STABILITY AND STORAGE

**Conditions:** protect from light (R2) close immediately after use

#### Substrate Start:

**Storage:** at 2 – 8°C

**Stability:** up to the expiration date

#### Sample Start (Working Reagent):

**Stability:** at 15 – 25°C 8 hours

at 2 – 8°C 5 days

protect from light!

Minimum allowable absorbance of the Working Reagent measured at 340 nm against water as reference is 1.1.

### SAMPLE STABILITY AND STORAGE

Loss of activity:	at 15 - 25°C	< 2% within 24 hours
	at 2 - 8°C	< 8 % within 3 days
Stability:	at -20 °C	at least 6 weeks

Discard contaminated specimens

### INTERFERING SUBSTANCES

No interference up to:

Ascorbic acid	30 mg/dl
Bilirubin	40 mg/dl
Triglycerides	2000 mg/dl

Hemoglobin interferes because LDH is released by erythrocytes.

### MANUAL TEST PROCEDURE

Bring reagents and samples to room temperature.

#### Substrate Start

Pipette into test tubes	25°C, 30°C	37°C
<b>Reagent 1</b>	1000 µl	1000 µl
<b>Sample</b>	20 µl	10 µl
Mix. Incubate for approximately 1- 5 min. Then add		
<b>Reagent 2</b>	250 µl	250 µl
Mix. Read initial absorbance against air after 1 minute and start a timer.		
Read absorbance again after exactly 1, 2 and 3 min. Determine ΔA/min. during the linear part of the assay.		

#### Sample Start

Pipette into test tubes	25°C, 30°C	37°C
<b>Working reagent</b>	1000 µl	1000 µl
<b>Sample</b>	20 µl	10 µl
Mix. Read initial absorbance against air after 1 minute and start a timer.		
Read absorbance again after exactly 1, 2 and 3 min. Determine ΔA/min. during the linear part of the assay.		

### CALCULATION (light path 1 cm)

LDH [U/L] = ΔA/min x Factor

**Factors:**

#### Substrate Start

	25°C or 30 °C	37°C
<b>Factor at 340 nm</b>	10080	20000
<b>Factor at 334 nm</b>	10275	20390
<b>Factor at 365 nm</b>	18675	37060

#### Sample Start

	25°C or 30 °C	37°C
<b>Factor at 340 nm</b>	8095	16030
<b>Factor at 334 nm</b>	8250	16345
<b>Factor at 365 nm</b>	15000	29705

### UNIT CONVERSION

U/L x 0.01667 = µkatal/L

### REFERENCE RANGE [U/L] \*

	25°C	30°C	37°C
<b>Adults</b>	< 240	< 346	< 480

\* It is recommended that each laboratory establishes its own normal range.

### TEST PRINCIPLE

Pyruvate + NADH + H<sup>+</sup> < LDH > Lactate + NAD<sup>+</sup>

Reaction is buffered at physiological pH to favor equilibrium to lactate.

### ABBREVIATIONS

LDH	= Lactate Dehydrogenase
NAD <sup>+</sup>	= Nicotinamide Adenine Dinucleotide
NADH	= Reduced NAD



## PERFORMANCE CHARACTERISTICS

### LINEARITY

The test has been developed to determine LDH activities which correspond to a maximal  $\Delta A/\text{min}$  of 0.15 at 340 and 334nm or 0.08 at 365nm.

If these values are exceeded the sample should be diluted 1 + 10 with NaCl (9 g/L sodium chloride in water) and results multiplied by 11.

### PRECISION (at 25 °C)

Intra-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	142	5.50	3.86
Sample 2	245	4.95	2.01
Sample 3	497	8.39	1.69

Inter-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	144	3.09	2.13
Sample 2	248	4.53	1.82
Sample 3	492	6.23	1.26

### METHOD COMPARISON

A comparison between Diagnostica LDH-P (y) and a commercially available test (x) using 78 samples gave following results:  $y = 1.03x + 2.13$  U/L;  $r = 0.999$ .

### QUALITY CONTROL

All control sera with LDH values determined by this method can be used. We recommend:

DIA030012	1 x 5 ml	DIACON N	Assayed Control Serum Normal
DIA030022	1 x 5 ml	DIACON P	Assayed Control Serum Abnormal

### CALIBRATION

The use of a LDH Calibrator is optional. We recommend:

DIA040012	1 x 3 ml	DIACAL	Assayed Multi Calibration Serum
-----------	----------	--------	---------------------------------

### AUTOMATION

Special adaptations for automated analyzers can be made on request.

### Warnings and Precautions

1. The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
2. Take the necessary precautions for the use of laboratory reagents.

### WASTE MANAGEMENT

Please refer to local legal requirements.



